



Laboratory Service Manual

Lab Dept: Chemistry

Test Name: TESTOSTERONE, TOTAL AND FREE

General Information

Lab Order Codes: FTES

Synonyms: N/A

CPT Codes: 84402 – Testosterone, free
84270 – Sex hormone binding globulin

84403 – Testosterone, total
Test Includes: Total testosterone (ng/mL), Free Testosterone (pg/mL), % free testosterone, and sex hormone binding globulin (SHBG)(nmol/L)

Logistics

Test Indications: Most circulating testosterone is bound to a specific carrier protein usually called sex hormone binding hormone (SHBG) or testosterone-binding globulin (TeBG). A small portion of total testosterone exists in the free or unbound state available for entry into the cells of target organs. Free testosterone is probably the physiologically active component of testosterone. The serum free testosterone assay was developed because the measurement total testosterone does not always give an accurate picture of the amount of testosterone available to target organ cells. A number of conditions or medications are known to increase or decrease the TeBG concentration, which may cause total testosterone concentration to change but without necessarily influencing the free testosterone concentration, many of the influences known to alter thyroxine-binding proteins also influence TeBG in a similar way. Thus, estrogens increase the concentration of TeBG and androgens tend to decrease it.

Lab Testing Sections: Chemistry - Sendouts

Referred to: Esoterix, Inc. (Test# 500290)

Phone Numbers:

Minneapolis: 612-813-6280

Saint Paul: 651-220-6550

Test Availability: Daily, 24 hours



Laboratory Service Manual

Turnaround Time: 3 – 5 days, test set up Monday - Friday

Special Instructions: N/A

Specimen

Specimen Type: Whole blood

Container: Red top (plain, no gel) tube

Draw Volume: 6.0 mL (Minimum: 3.0 mL) blood

Processed Volume: 2.0 mL (Minimum: 1.0 mL) serum

Note: Submission of the minimum volume does not permit repeat analysis.

Collection: Routine venipuncture

Special Processing: Lab Staff: Centrifuge specimen within 1 hour of collection. Remove serum aliquot into a screw-capped round bottom plastic vial. Store and ship frozen in a plastic vial. Forward promptly.

Patient Preparation: None

Sample Rejection: Specimen not separated within 1 hour of draw.

Interpretive

Reference Range:

Total Testosterone:



Laboratory Service Manual

Premature Infants:				
Age	Males		Females	
26 - 28 weeks, day 4	59 - 125 ng/dL		5 - 16 ng/dL	
31 - 35 weeks, day 4	37 - 198 ng/dL		5 - 22 ng/dL	
Infants: Full Term Newborns:				
Males*		Females**		
75 - 400 ng/dL		20 - 64 ng/dL		
Infants: 1-7 months				
*Males: Levels decrease rapidly the first week to 20 - 50 ng/dL, then increase to 60 - 400 ng/dL between 20 - 60 days. Levels then decline to prepubertal range by 7 months.				
**Females: Levels decrease during the first month to <10 mg/dL and remain there until puberty.				
Prepubertal Children:				
1 - 10 years		<3 - 10 ng/dL		
Puberty:				
Tanner Stage	Males ng/dL	Age (years)	Females ng/dL	Age (years)
1	<3 - 10	<9.8	< 3 - 10	<9.2
2	18 - 150	9.8 - 14.5	7 - 28	9.2 - 13.7
3	100 - 320	10.7 - 15.4	15 - 35	10.0 - 14.4
4	200 - 620	11.8 - 16.2	13 - 32	10.7 - 15.6
5	350 - 970	12.8 - 17.3	20 - 38	11.8 - 18.6
Adult:		Male		Female



Laboratory Service Manual

20 – 50 years:	350 – 1030 ng/dL	10 – 55 ng/dL
Free Testosterone:		
Males	% Free	pg/dL
Full Term Infants Males		
1 – 15 days	0.9 – 1.7	1.5 – 31
1 – 3 months	0.4 – 0.8	3.3 – 18
3 – 5 months	0.4 – 1.1	0.7 – 14
5 – 7 months	0.4 – 1.0	0.4 – 4.8
Prepubertal Children Males		
1 – 10 years	0.4 – 0.9	0.15 – 0.6
Puberty Males	Comprehensive values are currently unavailable	
Adults Males	1.5 – 3.2	52 - 280
Females	% Free	pg/dL
Full Term Infants Females		
1 – 15 days	0.8 – 1.5	0.5 – 2.5
1 – 3 months	0.4 – 1.1	0.1 – 1.3
3 – 5 months	0.5 – 1.0	0.3 – 1.1
5 – 7 months	0.5 – 0.8	0.2 – 0.6
Prepubertal Children Females		
1 – 10 years	0.4 – 0.9	0.15 – 0.6
Puberty Females	Comprehensive values are currently unavailable.	
Adults Females	0.8 – 1.4	1.1 – 6.3



Laboratory Service Manual

Sex Hormone Binding Globulin		
Age	Males (nmol/L)	Females (nmol/L)
1 month - 2 years:	60 – 252	60 – 252
2 - 8 years:	72 – 220	72 - 220
9 - 18 years:	16 – 100	36 - 125
>18 years:	20 – 60	40 - 120

Critical Values: N/A

Limitations: N/A

Methodology: Binding Capacity/HPLC Tandem Mass Spectrophotometer

Contraindications: N/A

References: Esoterix Inc. "Expected Value and S.I. Unit Conversion Table", 5th Edition

Esoterix Laboratories Web Page www.esoterix.com