



Laboratory Service Manual

Lab Dept: Anatomic Pathology

Test Name: ANGELMAN MICRODELETION FISH (15Q11-13)

General Information

Lab Order Codes: ANSY

Synonyms: N/A

CPT Codes: 88271 x2 – Molecular cytogenetics; DNA probe, each
88273 – Chromosome in situ hybridization, analyze 10-30 cells
88275 – Interphase in situ hybridization, analyze 100 – 300 cells

Test Includes: Angelman and Prader Willi are both at 15q11-13. Please indicate if both should be tested, or if Angelman only is desired.

Logistics

Test Indications: N/A

Lab Testing Sections: Anatomic Pathology - Sendouts

Referred to: Abbott Northwestern Hospital Laboratory – Cytogenetics (Ultra Panel Code: ANG)

Phone Numbers:

Minneapolis: 612-813-6280

Saint Paul: 651-220-6550

Test Availability: Daily until 7:00 pm

Turnaround Time: 2 – 21 days

Special Instructions: Request form must include date and time of specimen collection, attending physician, diagnosis (or reason for referral), clinical status, information regarding medication or transfusions, and specimen type. Angelman and Prader Willi are both at 15q11-13.

Specimen

Specimen Type: Whole blood



Laboratory Service Manual

Container:	Green top (Na heparin) tube
Draw Volume:	5 mL (Minimum: 1.0 mL)
Processed Volume:	Same as Draw Volume
Collection:	Routine venipuncture
Special Processing:	Lab Staff: Do Not centrifuge specimen. Store and send at room temperature in original collection container. Use ANW request form B. Call (612-863-4541) prior to sending specimen. Send to MOB 402 (Allina Cytogenetics). Note: Order both probes if Angelman and Prader are ordered. If the care provider orders FISH: 15q11-13, call and clarify the order since 4 microdeletions are located in this area.
Patient Preparation:	None
Sample Rejection:	Frozen or refrigerated samples, mislabeled or unlabeled specimens

Interpretive

Reference Range:	Interpretive report
Critical Values:	N/A
Methodology:	Determined by Laboratory Director based on clinical information.
References:	Abbott Northwestern Cytogenetics Lab May 2008