

Centricity™ Universal Viewer Zero Footprint Client DICOM® Viewer

Tools, Icons and Keyboard Shortcuts Job Aid











Intended Audience: Radiology and Cardiology users

















































Applicable Foundations: Centricity PACS (CPACS), Centricity PACS-IW, Enterprise Archive (EA) and Edison True PACS












Toolbar Tools/Icons

Icon	Name	Description
	Study Search	Display the Study Search screen. Note: If launching via an external system, this option may be disabled.
	Navigator	Toggle Navigator open/closed to view structured report, key image notes or studies in the ZFP viewer.
	Series Selector	Toggle Series Selector open/closed to access comparison exams for the selected patient.
	Layout	Change the layout of the active viewport. The toolbar shows the currently displayed layout in the active viewport.
	Sheet View Layout	Hover over to create a sheet view layout in the active viewport. Once the sheet view is created, the sheet view controls display in the toolbar. 
	Scroll	Scroll through images in a series.
	Zoom	Zoom in and out on an image. The drop-down menu includes: <ul style="list-style-type: none"> Original Resolution: Resize the image to its original captured image size. Zoom to Fit: Resize the image to its maximum size inside the viewport without cropping. 2x: Resize the image two times the original resolution. 3x: Resize the image three times the original resolution. 4x: Resize the image four times the original resolution.
	Pan	Activate the pan tool.

Centricity Universal Viewer Zero Footprint DICOM Client Viewer Tools, Icons and Keyboard Shortcuts Job Aid










Icon	Name	Description																																				
	Window/level	Activate the window/level tool.																																				
	Invert	Invert the grayscale of the image in the active viewport.																																				
	Flip Rotate	Options in the drop-down menu include:  Rotate an image to the right.  Flip an image vertically.  Flip an image horizontally.																																				
	Annotations and Measurements	Options in the drop-down menu include: <table border="1"> <tbody> <tr> <td></td> <td>Distance</td> <td>Measure distance.</td> </tr> <tr> <td></td> <td>Angle</td> <td>Draw a normal angle or Cobb angle.</td> </tr> <tr> <td></td> <td>Elliptical Region</td> <td>Draw an elliptical region around a region of interest (ROI).</td> </tr> <tr> <td></td> <td>Ellipse</td> <td>Draw an ellipse (oval) around a region of interest.</td> </tr> <tr> <td></td> <td>Point Measurement</td> <td>Measure the pixel value of a point (Hounsfield units).*</td> </tr> <tr> <td></td> <td>Freehand Region</td> <td>Draw a region of interest and the measurement value is displayed.**</td> </tr> <tr> <td></td> <td>Freehand</td> <td>Draw a freehand region of interest <u>without</u> a measurement value.*</td> </tr> <tr> <td></td> <td>Rectangle</td> <td>Draw a rectangle around a region of interest.</td> </tr> <tr> <td></td> <td>Text</td> <td>Create a text annotation.</td> </tr> <tr> <td></td> <td>Arrow</td> <td>Add an arrow pointer as an annotation.</td> </tr> <tr> <td></td> <td>CTR Measurement</td> <td>Perform a CTR measurement of the heart and chest.</td> </tr> <tr> <td></td> <td>Rectangle Shutter</td> <td>Neutralize (exclude) the display of any pixels located outside of a shutter shape.</td> </tr> </tbody> </table> <p>*This tool is only available in version 6.0 SP11.5 and higher. **For version 6.0 SP11.5, this tool was renamed to "Freehand Region". In versions prior to 6.0 SP11.5, this tool was called "Freehand". No changes were made to its functionality.</p>		Distance	Measure distance.		Angle	Draw a normal angle or Cobb angle.		Elliptical Region	Draw an elliptical region around a region of interest (ROI).		Ellipse	Draw an ellipse (oval) around a region of interest.		Point Measurement	Measure the pixel value of a point (Hounsfield units).*		Freehand Region	Draw a region of interest and the measurement value is displayed.**		Freehand	Draw a freehand region of interest <u>without</u> a measurement value.*		Rectangle	Draw a rectangle around a region of interest.		Text	Create a text annotation.		Arrow	Add an arrow pointer as an annotation.		CTR Measurement	Perform a CTR measurement of the heart and chest.		Rectangle Shutter	Neutralize (exclude) the display of any pixels located outside of a shutter shape.
	Distance	Measure distance.																																				
	Angle	Draw a normal angle or Cobb angle.																																				
	Elliptical Region	Draw an elliptical region around a region of interest (ROI).																																				
	Ellipse	Draw an ellipse (oval) around a region of interest.																																				
	Point Measurement	Measure the pixel value of a point (Hounsfield units).*																																				
	Freehand Region	Draw a region of interest and the measurement value is displayed.**																																				
	Freehand	Draw a freehand region of interest <u>without</u> a measurement value.*																																				
	Rectangle	Draw a rectangle around a region of interest.																																				
	Text	Create a text annotation.																																				
	Arrow	Add an arrow pointer as an annotation.																																				
	CTR Measurement	Perform a CTR measurement of the heart and chest.																																				
	Rectangle Shutter	Neutralize (exclude) the display of any pixels located outside of a shutter shape.																																				
	Overlay	Show or hide the image overlays.																																				
	Cross Reference	Turn on or turn off the cross-reference lines.																																				
	Presentation State	Apply a different presentation state or reset a presentation state.																																				

Centricity Universal Viewer Zero Footprint DICOM Client Viewer Tools, Icons and Keyboard Shortcuts Job Aid

Icon	Name	Description
	Bookmarked Images	Apply a bookmark from a drop-down list of bookmarks referenced in the current study. This icon is disabled when no bookmarks are referenced in a study.
	Synchronize	Synchronize parallel series. This icon is disabled for multi-frame images.
	External Application	Launch a web-based external application from the drop-down menu.
	Share Study URL	Copy or email a study URL. The email recipient can open the study by clicking the URL from the email. (Email recipient must be a valid ZFP user.)
	Reports/Notes	Display a report or clinical document associated with a study (if available for that study).
	Reset	Remove all annotations and manipulations from the active viewport.
	MIP/MPR	Hover over this icon to construct standard MPR, MIP, or volume rendering (VR) images for a study.
	Logout	Log out of ZFP.
	Lossy/Lossless Image Quality	<p>When this icon is selected, lossless progression is disabled for image display. Images are loaded as lossy only.</p> <p>When this icon is not selected, lossless progression is enabled for image display. This is the default.</p> <p>Note: The tooltip shows the current network bandwidth. ZFP tests your network speed after downloading all lossy images for a study. If your network speed is 25 Mbps or lower, then ZFP selects this icon while the images are loading so that you can view images faster.</p>
	About	Open the About dialog box. This icon is for the ZFP desktop version only.
	Help	<p>Options in the drop-down menu include:</p> <ul style="list-style-type: none"> ○ User Manual - Open the ZFP online help. ○ System Prerequisites - Display the System Prerequisites dialog box. This dialog box contains the optimal ZFP system prerequisites related to browsers, operating systems, display size, and product model. ○ Help Guide - Display static ZFP UI images with tips that highlight basic areas of functionality, such as how to navigate the Study Search screen, and tips for image navigation and manipulation. Not supported on mobile phones.

Icon	Name	Description
		<ul style="list-style-type: none"> ○ Security and Privacy - Display the Security and Privacy screen. This informational screen is intended to remind users about their privacy and security obligations

Keyboard Shortcuts

Keyboard Shortcut		Action
	Press and hold right arrow key	Scroll through the next images
	Press and hold left arrow key	Scroll through the previous images
	Press down arrow key	Next series
	Press up arrow key	Previous series
	Press Insert key	Middle image in series
	Press Home key	First image in series
	Press End key	Last image in series
	Press right arrow key	Next image in series
	Press left arrow key	Previous image in series

GE HealthCare reserves the right to make changes in specifications and features, or discontinue the product or service described at any time, without notice or obligation. This does not constitute a representation or warranty or documentation regarding the product featured. All illustrations or examples are provided for informational or reference purposes and/or as fictional examples only. Your product features and configuration may be different than those shown. Information contained herein is proprietary to GE HealthCare. No part of this publication may be copied or reproduced for any purpose without prior written permission from GE HealthCare.