Lab Dept: Anatomic Pathology

Test Name: FISH LOCUS SPECIFIC (METAPHASE)

**General Information** 

Lab Order Codes: FISUM

**Synonyms:** Fluorescent in Situ Hybridization; FISH Select Locus; FISH Constitutional;

FISH Congenital

88273 - Chromosomal in situ hybridization, analyze 10-30 cells

**CPT Codes:** 88271 – Molecular cytogenetics, DNA probe, each

88291 - Cytogenetics and molecular cytogenetics, interpretation and report

88230 - Tissue culture for non-neoplastic disorders; lymphocyte

Charges may vary with type of probe, analysis and interpretation required, as determined by the testing laboratory. Each probe ordered will accrue a

probe and analysis charge.

Logistics

**Test Indications:** Detection of common genetic aberrations in patients that are of diagnostic

and prognostic significance. Results of this test must always be interpreted in the context of relevant clinical and pathological data and should not be

used alone for diagnosis.

**Lab Testing Sections:** Anatomic Pathology – Sendouts

Referred to: MHealth Fairview University of Minnesota Medical Center Cytogenetics

UM Test Code: LAB4736 (aka LAB1402) or FISHCG

Phone Numbers: MIN Lab: 612-813-6280

STP Lab: 651-220-6550

**Test Availability:** Daily

Turnaround Time: Results in 29 days

**Special Instructions:** 

Special tube required. See Container, contact the laboratory for appropriate

tube.

Testing laboratory requires the following additional information in the order:

Family History

Clinical Indications for Testing

Probe requested

Specimen

Specimen Type: Whole blood

Container: Green top (Sodium Heparin), no gel, tube

**Draw Volume:** 5 mL (Minimum: 3 mL) blood

**Processed Volume:** Same as draw volume

**Collection:** Routine venipuncture

**Special Processing:** Lab Staff: Do not centrifuge. Forward specimen to Sendouts. Store and ship

at room temperature. For optimal testing results, the specimen must arrive at reference lab within 24 hours of collection. Specimens may be viable for

the cell culture portion for up to 72 hours. Partial results from DNA

extraction may be available for specimens up to 7 days.

Specimens must be received in the Cytogenetics Laboratory Mon-Fri by 5:30 pm, weekends and holidays by 4:30 pm. Specimens received after

these cutoffs will be processed the following day.

Sendouts must follow prompts or manually enter the answers to the

required order questions in SmarTerm prior to Atlas entry.

Patient Preparation: None

Sample Rejection: Improper storage, mislabeled or unlabeled specimen. Clotted or frozen

specimen. Incorrect specimen collection tube (anticoagulant).

Interpretive

**Reference Range:** An interpretive report will be provided

Critical Values: N/A

Limitations: N/A

Methodology: Fluorescent in Situ Hybridization (FISH)

Fairview Diagnostic Laboratories January 2023 References:

2/5/2013: CPT update **Updates:** 

5/15/2019: Updated locus list 1/18/2023: Updated test name, reference lab codes, cpt codes and required

ordering information