Lab Dept: Urine/Stool

Test Name: 2,3-DINOR -11 BETA PROSTAGLANDIN F2

ALPHA, RANDOM URINE

General Information

Lab Order Codes: 23BPG

Synonyms: Prostaglandin, Mastocytosis

CPT Codes: 84150 – Prostaglandin, each

Test Includes: Prostaglandin reported in pg/mg Cr and Creatinine reported in mg/dL.

Logistics

Test Indications: Useful for screening mast cell activation disorders including systemic

mastocytosis using random urine specimens.

Lab Testing Sections: Urine/Stool - Sendouts

Referred to: Mayo Clinic Laboratories (MML Test: 23BPR)

Phone Numbers: MIN Lab: 612-813-6280

STP Lab: 651-220-6550

Test Availability: Daily, 24 hours

Turnaround Time: 2 – 6 days

Special Instructions: N/A

Specimen

Specimen Type: Urine

Container: Random urine cup

Draw Volume: 5 mL (Minimum: 4 mL) urine

Processed Volume: 5 mL (Minimum: 4 mL) urine

Collection: Random urine

Special Processing: Lab Staff:

Random Urine: Measure urine aliquot and place in Mayo Urine

transport tube (Supply T465) with no preservatives preferred. Store and

ship specimen at refrigerated temperatures.

Patient Preparation: Patients taking aspirin or nonsteroidal anti-inflammatory drugs (NSAID)

may have decreased concentrations of prostaglandin F2 alpha. If possible, discontinue for 2 weeks or 72 hours, respectively, prior to

collecting a specimen.

Sample Rejection: Mislabeled or unlabeled

Interpretive

Reference Range: All ages: <1802 pg/mg creatine

Critical Values: N/A

Limitations: Patients taking aspirin or nonsteroidal anti-inflammatory drugs (NSAID)

may have decreased concentrations of prostaglandin F2 alpha. If possible, discontinue for 2 weeks or 72 hours, respectively, prior to

collecting a specimen.

Methodology: 23BPG: Liquid Chromatography-Tandem Mass Spectrometry (LC-

MS/MS)

Creatintine: Enzymatic Colorimetric Assay

References: Mayo Clinic Laboratories June 2021

Updates: 6/22/2021: Updated to Random Urine at Mayo.