
Lab Dept: Chemistry

Test Name: FREE FATTY ACIDS

General Information

Lab Order Codes: FFAT

Synonyms: Free Fatty Acids, Total, Serum

CPT Codes: 82725 – Fatty acids; nonesterified

Test Includes: Free fatty acids, total level reported in mmol/L.

Logistics

Test Indications: Evaluation of metabolic status of persons with endocrinopathies.
Detection of pheochromocytoma and of glucagon, thyrotropin, and adrenocorticotropin secreting tumors.
Monitoring of control of diabetes mellitus.

Lab Testing Sections: Chemistry - Sendouts

Referred to: Mayo Medical Laboratories (MML Test: NEFA)

Phone Numbers: MIN Lab: 612-813-6280

STP Lab: 651-220-6550

Test Availability: Daily, 24 hours

Turnaround Time: 1 – 3 days, test is set up Monday - Friday

Special Instructions: See [Patient Preparation](#)

Specimen

Specimen Type: Blood

Container: SST (Gold, marble or red) tube

Draw Volume: 3 mL (Minimum: 1.5 mL) blood

Processed Volume: 1 mL (Minimum: 0.5 mL) serum

| | |
|-----------------------------|---|
| Collection: | Routine venipuncture |
| Special Processing: | Lab Staff: Centrifuge specimen within 45 minutes of draw. Remove serum aliquot and freeze immediately in a screw-capped round bottom plastic vial. Store and ship frozen. Forward promptly. |
| Patient Preparation: | Patient must be on an overnight fast (preferably 12 – 14 hrs) Patient must not consume any alcohol for 24 hrs before the specimen is drawn. |
| Sample Rejection: | Warm specimens; gross hemolysis; gross lipemia; grossly icteric; mislabeled or unlabeled specimens |

Interpretive

| | | |
|-------------------------|-----------------|--------------------|
| Reference Range: | Age: | Range: |
| | 0-15 years | Not established |
| | > or = 16 years | 0.00 – 0.72 mmol/L |

Critical Values: N/A

Limitations: Patient should fast for 12 to 14 hours; however, in prolonged fasting or starvation, free fatty acid levels rise as much as threefold.

Patient should abstain from alcohol for at least 24 hours.

In order to eliminate the generation of free fatty acids from triglycerides by serum lipases (causing erroneous elevations), serum should be frozen soon after it is drawn and shipped frozen.

Methodology: Enzymatic Colorimetric

References: [Mayo Medical Laboratories Web Page](#) November 2017

Updates: 6/14/2010: Reference range update, previously listed as <730 µEq/L.
11/14/2017: Collection container update.