
Lab Dept: Chemistry

Test Name: VITAMIN C

General Information

Lab Order Codes: VITC

Synonyms: Ascorbic Acid

CPT Codes: 82180 – Ascorbic Acid (Vitamin C), blood

Test Includes: Vitamin C (Ascorbic acid) level reported in mg/dL.

Logistics

Test indications: Diagnosing Vitamin C deficiency.
Used as an aid to deter excessive intake.

Lab Testing Sections: Chemistry - Sendouts

Referred to: Mayo Clinic Laboratories (Mayo test: VITC)

Phone Numbers: MIN Lab: 612-813-62

STP Lab: 651-220-6550

Test Availability: Daily, 24 hours

Turnaround Time: Performed Monday – Friday, results in 3-5 days

Special Instructions: Specimen needs to be placed on wet ice after collecting and processed in the lab within 4 hours of collection.

Specimen

Specimen Type: Blood

Container: Green (Lithium heparin) tube

Draw Volume: 3 mL (Minimum: 1.5 mL) blood

Processed Volume: 1 mL (Minimum: 0.5 mL) plasma

Collection: Routine blood collection

Special Processing:	Lab Staff: Specimen should come to lab on wet ice and must be processed with 4 hours of collection. Centrifuge specimen at 4 degrees C, remove plasma aliquot and place in Amber Vial (Mayo supply T192) to PROTECT FROM LIGHT. Store and ship at frozen temperatures.
Patient Preparation:	Fasting overnight (12-14 hours) (infants-draw prior to next feeding). Water can be taken as needed.
Sample Rejection:	Gross hemolysis; improper specimen; mislabeled or unlabeled specimens

Interpretive

Reference Range: 0.4 – 2.0 mg/dL

- * Values below 0.2 mg/dL indicate significant deficiency.
- * Values greater than or equal to 0.2 mg/dL and less than 0.4 mg/dL are consistent with a moderate risk of deficiency due to inadequate tissue stores.
- * Values of 0.4 to 2.0 mg/dL indicate adequate supply.

The actual level at which vitamin C is excessive has not been defined. Values above 3.0 mg/dL are suggestive of excess intake. Whether vitamin C in excess is indeed toxic continues to be uncertain. However, limited observations suggest that this condition may induce uricosuria and, in individuals with glucose-6 phosphate dehydrogenase deficiency, may induce increased red blood cell fragility.

Critical Values: N/A

Limitations: Testing in non-fasting specimens or the use of vitamin supplementation can result in elevated plasma vitamin concentrations. Reference values were established in patients who were fasting.

Methodology: Liquid Chromatography-Tandem Mass Spectrometry (LC-MS/MS)

References: [Mayo Clinic Laboratories](#) (August 2020)