
Lab Dept: Microbiology/Virology

Test Name: CMV BY PCR, BLOOD

General Information

Lab Order Codes: CMVQU

Synonyms: Cytomegalovirus PCR, Plasma; CMV DNA Detect/Quant, Plasma

CPT Codes: 87497 - Infectious agent detection by nucleic acid (DNA or RNA); cytomegalovirus amplified probe technique

Test Includes: CMV detection by PCR. This test **does not** include CMV Culture.

Logistics

Test Indications: Early detection of CMV viremia. Monitoring CMV disease progression and response to antiviral therapy.

Lab Testing Sections: Microbiology/Serology - Sendouts

Referred to: Mayo Clinic Laboratories (Test: CMVQN)

Phone Numbers: MIN Lab: 612-813-6280

STP Lab: 651-220-6550

Test Availability: Daily, 24 hours

Turnaround Time: 1 – 5 days, test set-up Monday - Saturday

Special Instructions: N/A

Specimen

Specimen Type: Blood

Container: Lavender top (EDTA) tube

Draw Volume: 4.5 mL (Minimum: 1.8 mL) blood

Processed Volume: 1.5 mL (Minimum: 0.6 mL) plasma

Collection: Routine blood collection, invert gently to mix.

Special Processing: Lab Staff: Centrifuge specimen within 2 hours of collection, remove plasma aliquot into a screw-capped plastic vial. Store and ship at frozen temperatures. Forward promptly.

Plasma specimen is stable frozen (preferred) for 84 days and refrigerated for 6 days.

Patient Preparation: None

Sample Rejection: Specimens other than plasma; mislabeled specimens or unlabeled specimens

Interpretive

Reference Range: Undetected

The quantification range of this assay is 35 to 10,000,000 IU/mL (1.54 log to 7.00 log IU/mL), with a 95% or higher limit of detection at 35 IU/mL.

A result of "Undetected" indicates the absence of cytomegalovirus (CMV) DNA in the plasma (see Limitations).

A result of "<35 IU/mL (<1.54 log IU/mL)" indicates that CMV DNA is detected in the plasma, but the assay cannot accurately quantify the CMV DNA present below this level.

A quantitative value (reported in IU/mL and log IU/mL) indicates the level of CMV DNA (ie, viral load) present in the plasma.

A result of ">10,000,000 IU/mL (>7.00 log IU/mL)" indicates that CMV DNA level present in plasma is above 10,000,000 IU/mL (7.00 log IU/mL), and the assay cannot accurately quantify CMV DNA present above this level.

Limitations: Cytomegalovirus (CMV) viral load results generated with this assay may be higher (up to 1.00 log IU/mL) than those from the previous Cobas Ampliprep/Cobas Taqman CMV Test (Roche Molecular Systems, Inc.), due to the differences in sensitivity of the assays.

Mutations within the highly-conserved regions of the CMV DNA polymerase (UL54) gene covered by cobas CMV may affect primers or probe binding resulting in the under-quantitation of virus or failure to detect the presence of virus. The cobas CMV assay mitigates this risk through the use of redundant CMV target sequence amplification primers.

Methodology: Real-Time Reverse Transcription-Polymerase Chain Reaction (RT-PCR)

Additional Information:

CMV may be present in the blood in healthy seropositive individuals at very low levels. CMV serology may confirm current infection.

References:

[Mayo Clinic Laboratories Web](#) September 2023

Updates:

3/18/2004: Test moved from Eastern Virginia Medical School Department of Pathology to Mayo Medical Laboratories.

6/3/2013: Method change, previously listed as Real-Time PCR/DNA Probe Hybridization

8/19/2015: Draw volume change

8/3/2013: Draw volume change, updated PCR method.

9/21/2023: Updated turnaround time, added specimen stability, corrected links.