
Lab Dept: Urine/Stool

Test Name: KETONES, URINE

General Information

Lab Order Codes: UKE

Synonyms: Urine Ketones; Nitroprusside Reaction for Ketones

CPT Codes: 81003 – Urinalysis, by dipstick; automated, without microscopy

Test Includes: Screen for urine ketones.

Logistics

Test Indications: Useful for evaluation of ketonuria, detection of acidosis, ketoacidosis, fasting, starvation, high protein diets, diabetes mellitus, and stress-hormone excess.

Lab Testing Sections: Urinalysis

Phone Numbers: MIN Lab: 612-813-6280

STP Lab: 651-220-6550

Test Availability: Daily, 24 hours

Turnaround Time: 1 hour

Special Instructions: Submit only one (1) of the following: Catheterized specimen or clean-catch specimen.

Note: Indicate type of specimen (catheterized) on request form. Date and time of collection are required on request form for processing. Transport specimen to laboratory immediately following collection

Specimen

Specimen Type: Urine

Container: Urine cup

Draw Volume: Entire specimen collection (catheter or clean catch)

Processed Volume: Minimum: 0.5 mL urine

Collection: A specimen collected by catheterization is optimal; however, a clean-catch or mid-stream specimen is also acceptable. Random, voided specimens will be accepted, but are the least desirable and are not recommended if a urine culture is also being requested.

Special Processing: N/A

Patient Preparation: None

Sample Rejection: Less than 0.5 mL urine submitted; mislabeled or unlabeled specimen

Interpretive

Reference Range: Negative

Critical Values: N/A

Limitations: Specimens containing large amounts of ascorbic acid or levodopa metabolites, valproic acid, phenozopyridine (Pyridium®), PSP dye, phenylketones or phthalein compounds may cause false-positives.

Methodology: Multistix®, Nitroprusside reaction

References: Strasinger SK (1989) Urinalysis and Body Fluids, FA Davis Company
Ringsrud KM, et al (1995) Urinalysis and Body Fluids: A Color Text and Atlas, Mosby
Brunzel N (1994) Fundamentals of Urine and Body Fluid Analysis, WB Saunders

Updates: 2/23/2012: Acetest is no longer performed as a confirmation of urine ketones.