


If this is the first time logging in with this device, you must first complete the steps in the [Remote access first time set up](#) guide.

Accessing GE PACS Web

Note: Access to GE PACS requires specific permissions. Please submit a Children's Service Desk Account Request for Remote Access to GE PACS, or call Children's IT Service Desk.

1. Log into Secure Access at <https://www.childrensmn.org/signin>
2. Click the link to **GE PACS Web**. 
3. Click the link for the **msxml.msi** install and then click **Run**.

Status	Check	Value	Description/Solution
ok	JavaScript	is enabled	
ok	Browser type/version	IE 7.0	
ok	Script Engine	11.0	
ok	Userdata persistence	is enabled	
✗	ActiveX	-	The installation of the following controls failed: - MSXML (please download the msxml.msi installer to install MSXML manually.) - Pictorial Index - Viewport

There is an error, Centricity Enterprise Web cannot go on.

Please refer to Chapter "Installing the Centricity Web Client" of the Installation Guide for further guidance or contact your system administrator.

[Click here to try again.](#)

Do you want to run or save **msxml.msi** (5.04 MB) from gepacsweb.childrensmn.org?

4. Continue through the installation and when prompted, click **Install Now**.

MSXML 4.0 SP2 Parser and SDK Setup

Welcome to the MSXML 4.0 SP2 Parser and SDK Setup Wizard

The Setup Wizard will install MSXML 4.0 SP2 Parser and SDK on your computer. Click Next to continue or Cancel to exit the Setup Wizard.

MSXML 4.0 SP2 Parser and SDK Setup

Choose Setup Type

Choose the setup type that best suits your needs

Install Now
Install MSXML 4.0 at C:\Program Files (x86)\MSXML 4.0

Customize...
Allows users to choose which program features will be installed and where they will be installed. Recommended for advanced users.

5. Install ActiveX, and the AMIViewer and PictorialIndex add-on modules.

Status	Check	Value	Description/Solution
	JavaScript	is enabled	
	Browser type/version	IE 7.0	
	Script Engine	11.0	
	Userdata persistence	is enabled	
	ActiveX	-	The installation of the following controls failed: - Pictorial Index - Viewport

There is an error, Centricity Enterprise Web cannot go on.

Please refer to Chapter "Installing the Centricity Web Client" of the Installation Guide for further guidance or contact your system administrator.

[Click here to try again.](#)

Internet Explorer blocked this website from installing an ActiveX control. [What's the risk?](#)

This webpage wants to run the following add-on: 'AMIViewer ActiveX Control Module' from 'GE Medical Systems (unverified publisher)'.

[What's the risk?](#)

This webpage wants to run the following add-on: 'PictorialIndex ActiveX Control Module' from 'GE Medical Systems (unverified publisher)'.

[What's the risk?](#)

Click the Centricity **Try Again** link and sign into GE PACS with your Children's MN username /password.