Lab Dept:

Urine/Stool

Test Name: SUCCINYLACETONE, URINE **General Information** Lab Order Codes: SAUU Synonyms: Succinylacetone, Random Urine **CPT Codes:** 83291 - Organic acids; single, quantitative **Test Includes:** Succinylacetone, urine creatinine and interpretation Logistics Test Indications: Diagnose and monitor patients with hereditary tyrosinemia type I. This compound is found only in patients with this disease, not other tyrosine disorders. Lab Testing Sections: Urine/Stool - Sendouts **Referred to:** Mayo Medical Laboratories (MML Test: ZW242) forward to ARUP Laboratories **Phone Numbers:** MIN Lab: 612-813-6280 STP Lab: 651-220-6550 **Test Availability:** Daily, 24 hours **Turnaround Time:** Varies. Results in 1-2 weeks **Special Instructions:** Submit a random urine collection. No preservative. Refrigerate specimen during and after collection. Please submit a Patient History for Biochemical Genetics that includes: age, gender, diet, drug therapy, and family history. Specimen **Specimen Type:** Urine, random **Container:** Plastic leakproof container (No preservative). Urine GUARD® collection container is preferred for a timed urine sample.

Draw Volume: Submit entire urine collection

Processed Volume:	5 mL (Minimum: 1 mL) urine
Collection:	A random urine (first morning specimens are preferred) sample may be obtained by voiding into a urine cup and is often performed at the laboratory. Bring the refrigerated container to the lab with appropriate labeling.
Special Processing:	Lab Staff: Mix collection well before aliquot is taken. Aliquot 5 mL (Minimum: 1 mL) urine. Store and ship at frozen temperatures. A frozen specimen is critical.
Patient Preparation:	None
Sample Rejection:	Mislabeled or unlabeled specimens; ambient specimens
Interpretive	
Reference Range:	0.000-0.300 mmol/mol CRT An interpretive report will be provided.
Critical Values:	N/A
Limitations:	N/A
Methodology:	Liquid Chromatography/Tandem Mass Spectrometry
References:	Mayo Medical Laboratories August 2014
	ARUP Laboratories August 2014