

## SmartZone Notifications


SmartZone provides noncritical notifications in a manner that does not disrupt your workflow, allowing you to view patient-relevant items and act on them as needed. Alerts will display in the Banner Bar and remain a notification until the item is addressed.

SmartZone alerts include the following:

- Passive notifications that do not require immediate action.  
**Note:** If a notification requires immediate action, a link will be present.
- Patient-specific information i.e., BMI, advanced care plans, behavioral care plans, etc.
- Information or notifications present when you open the chart or refresh it.

To address a SmartZone notification:


1. The notification displays as a link if you can take action, i.e. place an order, update a care plan, etc.
2. Click the link to address the notification. The appropriate section of the chart displays.
3. Click collapse (>) to collapse the SmartZone panel.





## Notifications

There are four categories of notifications.

- **Information:** Indicates an information notification.



- **Warnings:** A warning indicates action is required. Instructions on how to complete these actions are included in the warning message. This may be a link or may direct you to another area of the patient's chart for review.



- **Alerts:** Brings attention to a result of concern.
- **Care Pathways:** Offers suggestions for care based on patient information.