

October 2014

## 3M PIV Securement Dressing

Approved use in outpatient procedures

The 3M peripheral IV (PIV) securement dressing has been evaluated and approved for use in Minnetonka and SPS outpatient departments in place of the Statlock for outpatient procedures only. However, should these patients need to be admitted to inpatient area you may see this product. Although this product looks the same as the current Tegaderm it has enhanced properties and is approved for use as an IV Securement device. In the event a patient is admitted to an inpatient unit, the 3M securement device may stay in place until an IV dressing change is needed as per Children's IV Insertion & Care clinical standard.

Specific features include:

- Reinforced border, deep notch, and tape strips work together to enlist the entire dressing in securement
- Adhesives hold strongly, manage moisture, and release easily
- Designed to maintain securement by preventing edge lift, flexing with patient movement, and managing moisture.



Further evaluation will be done within the inpatient setting in the near future to determine effectiveness for long term use. For questions regarding use of this product, contact Children's vascular access team: Denise Rucker, Judith Goebel or Dawn Rocheford or Vascular Access manager, Andrew Berndt.