



Add Custom Links to the Links Toolbar

You can add up to two custom links to your Links Toolbar. To create a custom link:


1. Navigate to the specific website to copy the URL for the custom link. Right-click on the website address and select **Copy**.



2. From the Links Toolbar in Cerner, click the Toolbar Options drop-down. Click **Add or Remove Buttons** and click **Customize**.




3. Click **Add**.



Customizing the Links Toolbar

4. Enter the display name for the link and paste the URL in the appropriate fields.
5. Click **OK**.



6. Click **Close**.
The custom link will be available on your Links Toolbar.

