## Lab Dept:

## **Anatomic Pathology**

## Test Name: BRONCHOALVEOLAR LAVAGE CELL COUNT AND DIFFERENTIAL

## **General Information**

Lab Order Codes:	Bronchoalveolar Lavage Cell Count
Synonyms:	BAL, cell count; Bronch washings, cell count; Bronchial brushing, cell count
CPT Codes:	89050 – Cell count 88108-26 – Smear and interpretation
Test Includes:	When ordered, nucleated cell count and manual differential, culture, special stains for microorganisms, other special stains.
Logistics	
Test Indications:	Useful for differential diagnosis and management of interstitial lung disease, sarcoid, histiocytosis, and hypersensitivity.
Lab Testing Section:	Anatomic Pathology and Hematology
Phone Numbers:	MIN Lab: 612-813-5875
	STP Lab: 651-220-6550
Test Availability:	Monday - Friday
Turnaround Time:	1 - 3 days
Special Instructions:	N/A
Specimen	
Specimen Type:	Bronchoalveolar washings
Container:	Screw capped, sterile vial
Draw Volume:	40 - 70 mL of bronchoalveolar washings
Processed Volume:	Varies depending on orders.
Collection:	Maintain sterility and forward promptly at ambient temperature only.

Special Processing:	Specify need for counts and stains on the Bronchoalveolar Lavage Form (available from lab)
Patient Preparation:	Physician preference
Sample Rejection:	Insufficient specimen; mislabeled or unlabeled specimen
Interpretive	
Reference Range:	Differential Macrophages: 80 - 90%
	Lymphocytes: 10%
	Nucleated Cell Count: 10 - 15 x 10 <sup>3</sup> /100 mL
Critical Values:	N/A
Limitations:	Wide variability in cell type and numbers are recovered. Standardized volumes and concentrations are not yet established.
Methodology:	Light microscopy
Contraindications:	Severe hypoxemia with impending respiratory failure.
References:	Rosai, Juan MD, Ackerman's Surgical Pathology, 8 <sup>th</sup> ed, vols 1 and 2, CV Mosby
	Sternberg, Steven (1997) Histology for Pathologists, 2 <sup>nd</sup> ed, Lippincott/Raven